


*H2 controls direct pattern recognition command of multifunction grips for past-generation bebionic\*, TASKA, ZEUS, Esper, MIA, or Psyonic Ability Hand where no powered wrist or powered elbow is present.*

*(\*with serial numbers BBHxxx00, hand will require digital unlock code available upon request from Ottobock)*

## PART 1: PHYSICAL CONNECTIONS




### NOTE

#### TASKA LOW PROFILE HAND

If using a Taska Low Profile, the Coapt H2 cable connects to the Taska Switch Block (supplied by Taska) as follows:

- Coapt cable connection marked '**Close**' must plug into Switch Block port marked '1'
- Coapt cable connection marked '**Open**' must plug into Switch Block port marked '2'
- Power/battery cable (not supplied by Coapt) must plug into Switch Block port marked '3'
- Coapt cable connection marked '**Digital**' must plug into two-pin pigtail coming out of Taska hand



## PART 2: SOFTWARE SETTINGS

### ELBOW

#### ESPIRE CLASSIC PLUS

For more information on connecting our Coapt System to this device please reference our Espire Classic Plus + Coapt Compatibilites Guide.

### HANDS

#### OTTOBOCK (old generation) BEBIONIC\*

Control Strategy: Dual Site  
Open / Close Strategy: Electrodes: A,B  
Control Options: Electrodes  
Control Response: Proportional

In ELECTRODE SCREEN:

- For both signal channels, set the low thresholds to about 10% of the numbered range and the high thresholds to about 90% of the numbered range

Remember to press "Send to Hand" after making any changes

#### PSYONIC ABILITY HAND

Ensure the "Control System" setting in the PSYONIC ABILITY HAND app is set to "COAPT"

Ensure all other hand settings in the PSYONIC ABILITY HAND app are set to the factory defaults

### HANDS (CONT.)

#### MIA HAND

Settings will automatically be set upon each power-on

#### TASKA

Settings will automatically be set upon each power-on. Use the QR code for more information in the Taska|Coapt Clinician Guide



#### ZEUS HAND

Input Option: Coapt

Speed Control Strategy: Proportional

#### ESPER HAND

Ensure the grip settings in the Esper app are set exactly as shown

