


# 11

## CLINICAL CONNECTIONS & SETTINGS GUIDE


*11 controls pronation/supination of powered wrists that have rotational controllers and open/close operation of prosthetic hands. Pattern recognition command of multifunction grips is not included.*

### PART 1: PHYSICAL CONNECTIONS

#### OTTOBOCK ELECTRIC WRIST ROTATOR (10S17) & MYOROTRONIC (13E205)



#### FILLAUER PROWRIST ROTATOR (5010056)



### PART 2: SOFTWARE SETTINGS

#### WRISTS

##### OTTOBOCK ELECTRIC WRIST ROTATOR WITH MYOROTRONIC

Ensure Ottobock MyoRotronic set to 'program 1' defaults.  
Refer to the MyoRotronic manual for instructions.

##### FILLAUER PROWRIST ROTATOR

###### WRIST SET UP

- Electric Wrist Enabled: Yes, Wrist Direction: A

###### RELAX DELAY

- Relax Timer: 0.1

###### SWITCH CONTROL

- Switch Rates A: 75, Switch Rates B: 75, Switch Window: 0.1, Fast Access: On, Co-Contraction: Off, Auto-Switch Back: No

#### HANDS

##### OTTOBOCK SENSORHAND, MYOHAND VARIPLUS, SYSTEM ELECTRIC GREIFER, OR SYSTEM ELECTRIC HAND

Ensure 'control mode'/'program' 1

##### ADAMS HAND

Ensure the "Control System" setting in the Adam's Hand is set to "Coapt"

#### HANDS (CONT.)

##### OTTOBOCK BEBIONIC

Control Strategy: Dual Site

Open / Close Strategy: Electrodes: A,B

Control Options: Electrodes

Control Response: Proportional

In ELECTRODE SCREEN:

- For both signal channels, set the low thresholds to about 10% of the numbered range and the high thresholds to about 90% of the numbered range

Remember to press "Send to Hand" after making any changes

##### FILLAUER PROHAND, PROETD/ETD2

###### HAND SET UP

- Input Type: EMG, Hand Filtering: Quick, Channel: Dual, Control: Differential, Hand Direction: A

###### Thresholds/Outputs

- A: 15, B: 15

###### Input Gains

- A: 5, B: 5

**OTHER COMPATIBLE HANDS SHOULD BE SET TO  
FACTORY DEFAULT FOR OPTIMAL OPERATION**